



Serving your health with Excellence

ELECTROSURGICAL INSTRUMENTS

High Quality Instruments



Index

European Bipolar Forceps	Page 01-07
European Irrigation Bipolar Forceps	Page 08
European Non-Stick Bipolar Forceps	Page 09-12
USA 2 Pin Non-Stick Bipolar Forceps	Page 13-14
USA 2 Pin Bipolar Forceps	Page 15-19
SilverGlide Non-Stick Bipolar Forceps	Page 20-23
Single use USA 2-Pin Bipolar Forceps	Page 24-26
Single use European Bipolar Forceps	Page 27
Monopolar Cables	Page 28
Bipolar Forceps Cables	Page 29
Diathermy Instruments	Page 30-31
Bipolar Artery Sealer	Page 32
Bipolar Electrodes	Page 33
Electrodes 2.4mm	Page 34
Conization Electrodes	Page 35
Electrodes 4.0mm	Page 36
Patient Return Plate	Page 37
Laparoscopic Electrodes	Page 38
Electro surgical Instruments for Gynecology	Page 39-42

European Bipolar Forceps

Bayonet
L= 22.0 cm
Tip 1.0 mm



PS-201

Bayonet (Ang.Up)
L= 16.0 cm
Tip 1.0 mm



PS-202

Bayonet (Ang.Up)
L= 18.0 cm
Tip 1.0 mm



PS-203

Bayonet (Ang.Up)
L= 20.0 cm
Tip 1.0 mm



PS-204

Bayonet (Ang.Down)
L= 20.0 cm
Tip 1.0 mm



PS-205

Bayonet (Ang.Up)
L= 18.0 cm
Tip 1.0 mm



PS-206

Bayonet (Ang.Up)
L= 20.0 cm
Tip 1.0 mm



PS-207

European Bipolar Forceps

Bayonet
L= 16.0 cm
Tip 1.0 mm



PS-208

Bayonet (Ang.Up)
L= 18.0 cm
Tip 1.0 mm



PS-209

Bayonet (Ang.Up)
L= 20.0 cm
Tip 1.0 mm



PS-210

Bayonet (Ang.Up)
L= 16.0 cm
Tip 1.0 mm



PS-211

Bayonet (Ang.Down)
L= 18.0 cm
Tip 1.0 mm



PS-212

Bayonet (Ang.Up)
L= 20.0 cm
Tip 1.0 mm



PS-213

McPherson (Str)
L= 9.0 cm
Tip 0.5 mm



PS-214

European Bipolar Forceps

Coaptation Nadler
L= 10.2 cm
Tip 0.5 mm



PS-215

Iris (Str.)
L= 10.2 cm
Tip 0.5 mm




PS-216

Iris (Str.)
L= 10.8 cm
Tip 1.0 mm



PS-217

Straight
L= 16.0 cm
Tip 1.0mm



PS-218

Straight
L= 18.0 cm
Tip 1.0 mm



PS-219

Straight
L= 20.0 cm
Tip 1.0 mm



PS-220


Bayonet
L= 22.0 cm
Tip 1.0 mm



PS-221


European Bipolar Forceps

Straight
L= 20.0 cm
Tip 1.0 mm



PS-222

Straight
L= 22.0 cm
Tip 1.0 mm



PS-223

Straight (Ang.)
L= 16.0 cm
Tip 1.0 mm



PS-224

Straight (Ang.)
L= 18.0 cm
Tip 1.0 mm



PS-225

Straight (Ang.)
L= 20.0 cm
Tip 1.0 mm



PS-226

Bayonet
L= 16.0 cm
Tip 1.0 mm



PS-227

Bayonet
L= 18.0 cm
Tip 1.0 mm



PS-228

European Bipolar Forceps

Straight (Ang.)
L= 16.0 cm
Tip 1.0 mm



PS-229

Straight (Ang.)
L= 18.0 cm
Tip 1.0 mm



PS-230

Straight (Ang.)
L= 20.0 cm
Tip 1.0 mm



PS-231

Straight (Ang.)
L= 22.0 cm
Tip 1.0 mm



PS-232

Bayonet
L= 16.0 cm
Tip 1.0 mm



PS-233

Bayonet
L= 18.0 cm
Tip 1.0 mm



PS-234

Bayonet
L= 20.0 cm
Tip 1.0 mm



PS-235

European Bipolar Forceps

McPherson (Crv.)
L= 11.5 cm
Tip 0.5 mm



PS-236

Jeweler # 5
L= 11.5 cm
Tip 0.5 mm



PS-237

Jeweler # 5
L= 11.5 cm
Tip 0.5 mm



PS-238

Jeweler # 5
L= 11.5 cm
Tip 0.5 mm



PS-239

Adson
L= 12.1 cm
Tip 1.0 mm



PS-240

Straight
L= 16.0 cm
Tip 1.0 mm



PS-241

Bayonet
L= 18.0 cm
Tip 1.0 mm



PS-242

European Bipolar Forceps

McPherson (Str.)
L= 9.0 cm
Tip 0.5 mm



PS-243

McPherson (Crv.)
L= 11.5 cm
Tip 0.5 mm



PS-244

Jeweler # 5
L= 11.5 cm
Tip 0.0 mm



PS-245

Adson
L= 12.1 cm
Tip 0.1 mm



PS-246

Jeweler # 4
L= 11.5 cm
Tip 0.5 mm



PS-247

Jeweler # 7
L= 11.5 cm
Tip 0.5 mm



PS-248

Coaptation
L= 10.8 cm
Tip 0.5 mm



PS-249

European Irrigation Bipolar Forceps

Straight
L= 16.0 cm
Tip 1.0mm



PS-250

Straight
L= 20.0 cm
Tip 1.0mm



PS-251

Bayonet
L= 16.0 cm
Tip 1.0 mm



PS-252

Bayonet
L= 20.0 cm
Tip 1.0 mm



PS-253

European Non-Stick Bipolar Forceps

McPherson (Str.)
L= 9.0 cm
Tip 0.5 mm



BS-254

McPherson (Crv.)
L= 11.5 cm
Tip 0.5 mm



BS-255

Jeweler # 5
L= 11.5 cm
Tip 0.5 mm



BS-256

Adson
L= 12.1 cm
Tip 0.1 mm



BS-257

Jeweler # 4
L= 11.5 cm
Tip 0.5 mm



BS-258

Jeweler # 7
L= 11.5 cm
Tip 0.5 mm



BS-259

Coaptation
L= 10.8 cm
Tip 0.5 mm



BS-260

European Non-Stick Bipolar Forceps

Coaptation Nadler

L= 10.2 cm

Tip 0.5 mm



BS-261

Iris (Str.)

L= 10.8 cm

Tip 0.5 mm



BS-262

Iris (Str.)

L= 10.8 cm

Tip 1.0 mm



BS-263

Straight

L= 18.0 cm

Tip 1.0mm



BS-264

Straight

L= 20.0 cm

Tip 1.0 mm



BS-265

Straight

L= 20.0 cm

Tip 1.0 mm



BS-266

Bayonet

L= 22.0 cm

Tip 1.0 mm



BS-267

European Non-Stick Bipolar Forceps

Bayonet (Ang.Up)
L= 16.0 cm
Tip 1.0 mm



BS-268

Bayonet (Ang.Up)
L= 18.0 cm
Tip 1.0 mm



BS-269

Bayonet (Ang.Up)
L= 20.0 cm
Tip 1.0 mm



BS-270

Bayonet (Ang.Down)
L= 16.0 cm
Tip 1.0 mm



BS-271

Bayonet (Ang.Up)
L= 18.0 cm
Tip 1.0 mm



BS-272

Bayonet (Ang.Up)
L= 20.0 cm
Tip 1.0 mm



BS-273

European Non-Stick Bipolar Forceps

Straight (Ang.)
L= 16.0 cm
Tip 1.0 mm



BS-274

Straight (Ang.)
L= 18.5 cm
Tip 1.0 mm



BS-275

Straight (Ang.)
L= 20.0 cm
Tip 1.0 mm



BS-276

Straight (Ang.)
L= 22.0 cm
Tip 1.0 mm



BS-277

Bayonet
L= 16.0 cm
Tip 1.0 mm



BS-278

Bayonet
L= 18.0 cm
Tip 1.0 mm



BS-279

Bayonet
L= 20.0 cm
Tip 1.0 mm



BS-280

USA 2 Pin Non-Stick Bipolar Forceps

Cushing
L= 18.0 cm
Tip 1.0 mm



BBS-313

Gerald Bayo
L= 17.8 cm
Tip 1.0 mm



BBS-314

Jansen Bayo
L= 18.0 cm
Tip 1.0 mm



BBS-315

Scollive Gre
L= 19.1 cm
Tip 1.5 mm



BBS-316

Yasargil
L= 19.8 cm
Tip 1.0 mm



BBS-317

Hardy Bayo
L= 21.5 cm
Tip 1.0 mm



BBS-318

USA 2 Pin Non-Stick Bipolar Forceps

McPherson (Str.)
L= 9.0 cm
Tip 0.5 mm



BBS-306

McPherson (Crv.)
L= 9.0 cm
Tip 0.5 mm



BBS-307

Jeweler # 4
L= 11.5 cm
Tip 0.5 mm



BBS-308

Jeweler # 4
L= 11.5 cm
Tip 0.5 mm



BBS-309

Adson
L= 12.1 cm
Tip 1.0 mm



BBS-310

Semkin
L= 14.0 cm
Tip 1.0 mm



BBS-311

Cushing
L= 16.0 cm
Tip 1.0 mm



BBS-312

USA 2 Pin Bipolar Forceps

Scollive Green-Wood
L= 19.1 cm
Tip 1.5 mm



PB-302

Yasargil
L= 19.8 cm
Tip 1.5 mm



PB-303

Hardy Bayonet
L= 12.5 cm
Tip 1.0 mm



PB-304

Hardy Bayonet (Ang.)
L= 21.5 cm
Tip 1.0 mm



PB-305

USA 2 Pin Bipolar Forceps

McPherson (Str.)
L= 9.0 cm
Tip 0.5 mm



PB-295

McPherson (Crv.)
L= 9.0 cm
Tip 0.5 mm



PB-296

Coaptation Nadler
L= 10.8 cm
Tip 0.5 mm



PB-297

Coaptation
L= 10.8 cm
Tip 0.5 mm



PB-298

Jeweler # 5
L= 11.0 cm
Tip 0.5 mm



PB-299

Jeweler # 4
L= 11.5 cm
Tip 0.5 mm



PB-300

Tenzel
L= 11.5 cm
Tip 0.5 mm



PB-301

USA 2 Pin Bipolar Forceps

Jeweler # 7
L= 11.5 cm
Tip 0.5 mm



PB-288

Iris (Str.)
L= 10.8 cm
Tip 0.5 mm



PB-289

Iris (Crv.)
L= 10.8 cm
Tip 1.0 mm



PB-290

Adson
L= 12.1 cm
Tip 1.0mm



PB-291

Semkin
L= 14.0 cm
Tip 1.0 mm



PB-292

Cushing
L= 16.0 cm
Tip 1.0 mm



PB-293

Cushing
L= 11.5 cm
Tip 1.0 mm



PB-294

USA 2-Pin Bipolar Forceps

Cushing (Gursh)
L= 18.0 cm
Tip 1.0 mm



PB-281

Gerald (Str)
L= 17.8 cm
Tip 1.0 mm



PB-282

Gerald (Str)
L= 17.8 cm
Tip 1.0 mm



PB-283

Gerald Bayonet
L= 17.8 cm
Tip 1.0 mm



PB-284

Jansen Bayonet
L= 15.0 cm
Tip 1.0 mm



PB-285

Jansen Bayonet
L= 18.0 cm
Tip 1.0 mm



PB-286

Malis
L= 17.8 cm
Tip 1.0 mm



PB-287

USA 2 Pin Bipolar Forceps

Semkin
Irrigation



PB-319

Gerald Bayonet
Irrigation



PB-320

Jansen Bayonet
Irrigation



PB-321

Scollive Green-Wood
Irrigation



PB-322

Hardy Bayonet
Irrigation



PB-323

SilverGlide Non-Stick Bipolar Forceps

Bayonet Fu
L= 18.0 cr
Tip 1.0 mn



RN-324

Bayonet Fu
L= 20.0 cr
Tip 1.0 mn



RN-325

Bayonet Fu
L= 22.0 cr
Tip 1.0 mn



RN-326

Bayonet (lr
L= 20.0 cr
Tip 1.0 mn



RN-327

Bayonet (lr
L= 22.0 cr
Tip 1.0 mn



RN-328

Bayonet Fu
L= 22.0 cr
Tip 1.0 mn



RN-329

Bayonet Fu
L= 22.0 cr
Tip 0.5 mn



RN-330

Bayonet Fu
L= 24.0 cr
Tip 1.0 mn



RN-331

SilverGlide Non-Stick Bipolar Forceps

McPherson (Str.)
L= 9.0 cm
Tip 0.5 mm



RN-332

McPherson (Crv.)
L= 9.0 cm
Tip 0.5 mm



RN-333

Jeweler # 5
L= 11.5 cm
Tip 0.5 mm



RN-334

Adson
L= 12.1 cm
Tip 1.0 mm



RN-335

Cushing
L= 16.0 cm
Tip 1.0 mm



RN-336

Jansen Bayonet
L= 18.0 cm
Tip 1.0 mm



RN-337

Jansen Bayonet
L= 20.0 cm
Tip 1.0 mm



RN-338

Bayonet Full Non-
L= 16.0 cm
Tip 1.0 mm



RN-339

SilverGlide Non-Stick Bipolar Forceps

Jeweler # 5
L= 11.5 cm
Tip 0.5 mm



RN-340

Jeweler # 5
L= 11.5 cm
Tip 0.5 mm



RN-341

Straight Full Non-Stick
L= 12.0 cm
Tip 1.0 mm



RN-342

Straight (Ang.)
L= 16.0 cm
Tip 1.0 mm



RN-343

Bayonet Full Non-Stick
L= 18.0 cm
Tip 1.0 mm



RN-344

Bayonet Full Non-Stick
L= 20.0 cm
Tip 1.5 mm



RN-345

Bayonet Full Non-Stick
L= 22.0 cm
Tip 1.5 mm



RN-346

Bayonet Full Non-Stick
L= 22.0 cm
Tip 1.0 mm



RN-347

SilverGlide Non-Stick Bipolar Forceps

Jeweler # 5
L= 11.5 cm
Tip 0.5 mm



RN-355

Jeweler # 5
L= 11.5 cm
Tip 0.5 mm



RN-356

Straight Full Non-Stick
L= 12.0 cm
Tip 1.0 mm



RN-357

Straight (Ang.)
L= 16.0 cm
Tip 1.0 mm



RN-358

Bayonet Full Non-Stick
L= 18.0 cm
Tip 1.0 mm



RN-359

Bayonet Full Non-Stick
L= 20.0 cm
Tip 1.5 mm



RN-360

Bayonet Full Non-Stick
L= 22.0 cm
Tip 1.5 mm



RN-361

Bayonet Full Non-Stick
L= 22.0 cm
Tip 1.0 mm



RN-362

Single use USA 2-Pin Bipolar Forceps

Gerald Straight
Single Use



YW-348

Jansen Bayonet
Irrigation



YW-349

Jansen Bayonet
Single Use



YW-350

Jeweler # 5
Single Use



YW-351

Adson Forceps
Single Use



YW-352

Cushing Forceps
Size: 16.0 cm
Tip: 1.0 mm



YW-353

Jansen Bayonet ang-up
Size: 16.0 cm
Tip: 1.0 mm



YW-354

Single use USA 2-Pin Bipolar Forceps

Jansen Bayonet
Size: 18.5 cm
Tip 1.0 mm



YW-363

Jeweler # 4
Size: 11.5 cm
Tip 0.5 mm



YW-364

Adson Forceps
Size: 12.1 cm
Tip 1.0 mm



YW-365

Jansen Bayonet
Size: 18.0 cm
Tip 1.0 mm



YW-366

Jansen Bayonet
Size: 20.0 cm
Tip 1.0 mm



YW-367

Single use USA 2-Pin Bipolar Forceps

Jeweler # 5
Single Use



BY-373

Straight Forceps
Single Use



BY-374

Bayonet Forceps
Single Use



BY-375

Straight Forceps
Single Use



BY-376

Straight Forceps (Ang.)
Single Use



BY-377

Single use European Bipolar Forceps

Bayonet Forceps
Size: 18.0 cm
Tip 1.0 mm



YW-368

Bayonet Forceps
Size: 20.0 cm
Tip 1.0 mm



YW-369

Straging Forceps
Size: 20.0 cm
Tip 1.0 mm



YW-370

Straight Forceps (Ang.)
Size: 14.0 cm
Tip 1.0 mm



YW-371

Straight Forceps
Size: 14.0 cm
Tip 1.0 mm



YW-372

Monopolar Cables



Silicone Monopolar Cables 10" (3.0m)
SB-1501
 Instrument Connector: Screw Type
 Generator Connector: Standard Bovie 8.0mm Plug
 Compatibility: Valleylab, Conmed, Bovie



Silicone Monopolar Cables 10" (3.0m)
SB-1502
 Instrument Connector: 4.0mm Male Fitting
 Generator Connector: Standard 4.0mm Banana Plug
 Compatibility: Valleylab, Conmed & All ES Generators



Silicone Monopolar Endoscopic Cable 10" (3.0m)
SB-1503
 Instrument Connector: 4.0mm Female Plug
 Generator Connector: Standard Bovie 8.0mm Plug
 Compatibility: Valleylab, Conmed, Bovie



Silicone Monopolar Endoscopic Cable 10" (3.0m)
SB-1504
 Instrument Connector: 4.0mm Female Plug
 Generator Connector: Standard 4.0mm Banana Plug
 Compatibility: Valleylab, Conmed & All ES Generators



Silicone Diathermy Cable 10" (3.0m)
SB-1505
 Instrument Connector: 4.8mm Male Fitting
 Generator Connector: Standard Bovie 8.0mm Plug
 Compatibility: Valleylab, Conmed, Bovie



Silicone Diathermy Cable 10" (3.0m)
SB-1506
 Instrument Connector: 4.8mm Male Fitting
 Generator Connector: Standard 4.0mm Banana Plug
 Compatibility: Valleylab, Conmed & All ES Generators

Bipolar Forceps Cables



Silicone Bipolar Cables 10" (3.0m)
SB 2501

Instrument Connector: USA 2-Pin Plug
Generator Connector: 4.0mm Banana Plug
Compatibility: Valleylab, Conmed, Codman, Olympus and Mbst Bipolar Generator



Silicone Bipolar Cables 10" (3.0m)
SB 2502

Instrument Connector: European Flat Plug
Generator Connector: Martin Type
Compatibility: Aesculap, Martin, Wolf, Berchold



Silicone Bipolar Cables 10" (3.0m)
SB 2503

Instrument Connector: European Flat Plug
Generator Connector: 4.0mm Banana Plug
Compatibility: Valleylab, Conmed, Codman, Olympus & Mbst Bipolar Generator



Silicone Bipolar Cables 10" (3.0m)
SB 2504

Instrument Connector: USA 2-Pin Plug
Generator Connector: Martin Type
Compatibility: Aesculap, Martin, Wolf, Berchold



Silicone Bipolar Cables 10" (3.0m)
SB 2505

Instrument Connector: USA 2-Pin Plug
Generator Connector: Erbe Type
Compatibility: Eibe, Siemens, Downs



Silicone Bipolar Cables 10" (3.0m)
SB 2506

Instrument Connector: USA 2-Pin Plug
Generator Connector: Pin Distance 29.0mm
Compatibility: Valleylab, EMS, Dolley

Silicone Bipolar Cables 10" (3.0m)
SB 2507

Instrument Connector: European Flat Plug
Generator Connector: Pin Distance 29.0mm
Compatibility: Valleylab, EMS, Dolley



Diathermy Instruments



**Turnur Warwick Scissor
(Straight)**
SB-05-001



**Turnur Warwick Scissor
(Curved)**
SB-05-002



**Surgical Scissor
(Straight)**
SB-05-011



**Surgical Scissor
(Curved)**
SB-05-012



Abbey Sub Needle
17.5cm
SB-05-015



**Willson Hey Forceps
(Straight)**
SB-05-020

**Richie Forceps
(Straight)**
SB-05-038



**Willson Hey Forceps
(Curved)**
SB-05-022



**Artery Forceps
(Straight)**
SB-05-026



**Artery Forceps
(Curved)**
SB-05-028



31

Diathermy Instruments

Adson



PS-378

Mcldone



PS-379

Gillie



PS-380

Dissecting Forceps



Gerald (1x2)



PS-381

Adson Bayonet



PS-382

Adson Bayonet
Angle Up



PS-383

Remington Hob



Bipolar Artery Sealer



Artery Forceps Bi-Clamp

SB-18-004-25



Artery Forceps Non Stick

SB-18-006-25




Artery Forceps

SB-17-004-25



Artery Forceps

SB-17-001-25



Silicone Bi-Clamp Connecting cord

Bipolar Electrodes



Bi-Turbinate Probe
Size: 17.5 cm
SB-14-028

Electrodes 2.4mm



Blade Electrode
Disposable: 03-001



Ball Electrode (3.0 mm)
Disposable: 03-007



Circle Loop Electrode 4.0 mm
Reusable: 04-001



Ball Electrode (crv.) (3.0 mm)
Disposable: 03-008



Blade Electrode (crv.)
Disposable: 03-002



Circle Loop Electrode 4.0 mm
Reusable: 04-001



Needle Electrode
Disposable: 03-003



Ball Electrode (4.0 mm)
Disposable: 03-009



Circle Loop Electrode 8.0 mm
Reusable: 04-003



Needle Electrode (crv.)
Disposable: 03-004



Ball Electrode (crv.) (4.0 mm)
Disposable: 03-010



Diamond Electrode 8.0 mm
Reusable: 04-002



Ball Electrode (2.0 mm)
Disposable: 03-005



Ball Electrode (5.0 mm)
Disposable: 04-011



Ball Electrode (crv.) (2.0 mm)
Disposable: 03-006




Ball Electrode (crv.) (5.0 mm)
Disposable: 03-012



Tip Cleaner 03-35



Blade Electrode (Shaft 13.0 cm)
Disposable: 03-013



Holster 038



Blade Electrode (Shaft 13.0 cm)
Disposable: 03-015



Blade Electrode (Shaft 13.0 cm)
Disposable: 03-023



Blade Electrode (Shaft 13.0 cm)
Disposable: 03-007




Conization Electrode
Disposable: 03-025


Conization Electrodes




Conization Electrode
Disposable #03-026




Conization Electrode
Disposable #03-027



Conization Electrode
Disposable #03-029



Conization Electrode
Disposable #03-030



Conization Electrode
Disposable #03-030



Extender Electrode
03-032 Size: 7.0 cm



Extender Electrode
#03-033 Size:13.0 cm



Extender Electrode
#03-034 Size:18.0 cm

Electrodes 4.0mm



Blade Electrode
Reusable 20-001



Blade Electrode (Crv)
Reusable 20-002



Needle Electrode
Reusable 20-003



Needle Electrode (Crv)
Reusable 20-004



Ball Electrode (2mm)
Reusable 20-005



Ball Electrode (Crv) (2mm)
Reusable 20-006



Blade Electrode (Shaft 13.0cm)
Reusable 20-018



Needle Electrode (Shaft 13.0cm)
Reusable 20-019



Ball Electrode 3mm (Shaft 13.0 cm)
Reusable 20-020




Ball Electrode 4mm (Shaft 13.0cm)
Reusable 20-021



Ball Electrode 5mm (Shaft 13.0cm)
Reusable 20-022



Ball Electrode
Reusable 20-007



Ball Electrode (Crv) 3mm
Reusable 20-008



Ball Electrode (4mm)
Reusable 20-009



Ball Electrode (Crv) (4mm)
Reusable 20-010



Ball Electrode (5mm)
Reusable 20-011



Ball Electrode (Crv) (5mm)
Reusable 20-012



Round Loop Electrode
(7.0cm) 4.0mm
Reusable 20-013



Round Loop Electrode
(7.0cm) 6.0mm
Reusable 20-014



Round Loop Electrode
(7.0cm) 8.0mm
Reusable 20-015



Round Loop Electrode
(7.0cm) 10.0mm
Reusable 20-016




Plate Electrode Small
Reusable 20-017





Plate Electrode Large
Reusable 20-018

37

Patient Return Plate



Patient Return Plate Reusable
5m Cable #25-001

Patient Return Plate Reusable
(5m Cable) 25-002

Patient Return Cable
25-006-a (5m cable)
25-006-b (3m Cable)

Patient Return Cable
25-007-a (5m cable)
25-007-b (3m Cable)

S.S Patient Return Plate
(16x12cm)25-003

Laparoscopic Electrodes

WITHOUT SUCTION



15-001 L-Hook Electrode (36cm)



15-002 J-Hook Electrode (36cm)



15-003 Needle-Hook Electrode (36cm)



15-004 Spatula-Hook Electrode (36cm)



15-005 Ball-Hook Electrode (36cm)

WITH SUCTION



15-006 L-Hook Electrode (Irrigation) (36cm)



15-007 J-Hook Electrode (Irrigation) (36cm)



15-008 Needle-Hook Electrode (Irrigation) (36cm)



15-009 Spatula-Hook Electrode (Irrigation) (36cm)



15-010 Ball-Hook Electrode (Irrigation) (36cm)

Laparoscopic Electrodes



Electro surgical Instruments for Gynecology



Lateral Vaginal
Retractor - #22-001



Kogan Endocervical
Speculum - #22-003



Lateral Vaginal Retractor
(wide-Open) #22-006



Rigby Retractor
22-007



Rigby Retractor
22-005



3mm 130-133



4mm 130-134



3mm 130-133

Electro surgical Instruments for Gynecology




Weisman Speculum with Suction (Left Open)
Small: 21-061
Medium: 21-062
Large: 21-063



Weisman Speculum with Suction (Right Open)
Small: 21-071
Medium: 21-072
Large: 21-073



Dressing Forceps
23-001 (21.0 cm)
23-001 (25.0 cm)



Tissue Forceps
23-003 (21.0 cm)
23-004 (25.0 cm)



Iris Hook
23-005 (25.0 cm)



Emmet Hook
23-006 (25.0 cm)



Straight Hook
23-007 (25.0 cm)



Two-Prong Hook
23-008 (25.0 cm)

Electro surgical Instruments for Gynecology



Graves Vaginal Speculum
with smoke evacuation tube
Small: 21-001
Medium: 21-002
Large: 21-003



VU-Max Speculum
with smoke evacuation tube
Small: 21-011
Medium: 21-012
Large: 21-013



Pederson Vaginal Speculum
with smoke evacuation tube
Small: 21-041
Medium: 21-042
Large: 21-043



VU-More (Graves) Speculum
with smoke evacuation tube
Small: 21-031
Medium: 21-032
Large: 21-033

Cusco Speculum
with smoke evacuation tube
Small: 21-021
Medium: 21-022
Large: 21-023



Collin Speculum
with smoke evacuation tube
Small: 21-051
Medium: 21-052
Large: 21-053



Electro surgical Instruments for Gynecology




**Kogan Endocervical
Speculum**
SB-#25-003



**Lateral Vaginal
Retractor**
SB-#25-001



**Graves Speculum
with smoking evacuation tube**
Small: 21-001
Medium: 21-002
Large: 21-003



**Pederson Graves Speculum
with smoking evacuation tube**
Small: 21-001
Medium: 21-002
Large: 21-003



**Collin Speculum
with smoking evacuation tube**
Small: 24-051
Medium: 24-052
Large: 24-053



SUBLIME SURGICAL

45 Waters Edge,
Canterbury CT1 1WX
United Kingdom

Tel: +447427677804

Web: www.sublimesurgical.com

E-mail: info@sublimesurgical.com