

PLASTIC SURGERY INSTRUMENTS

High Quality Instruments



INDEX

BREAST DISSECTORS -
ELEVATORS

Pg 1-5

BREAST MARKERS AND
MEASURING CALIPERS

Pg 6-8

BREAST RETRACTORS

Pg 9-10

ENDOSCOPIC FACE AND FOREHEAD
LIFT INSTRUMENTS

Pg 11-13

FACE LIFT SCISSORS

Pg 14-15

**AGRIS-DINGMAN Breast Dissectors
Left and Right Dissectors 36 cm**

SKU: SSPS 311

Instrument Type: Dissector
Length: 36cm
Working Surface Style Left and Right
Material: Japanese Stainless Steel
Reusable
Non-Sterile



**AGRIS-DINGMAN Breast Dissectors
Graduated Right 36 cm**

SKU: SSPS 312

Instrument Type: Dissector
Length: 36cm
Working Surface Style Right
Material: Japanese Stainless Steel
Reusable
Non-Sterile



**AGRIS-DINGMAN Breast Dissectors
Graduated Left 36 cm**

SKU: SSPS 313

Instrument Type: Dissector
Length: 36cm
Working Surface Style Left
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Breast Dissector Oval Spatulated Blade

SKU: SSPS 314

Instrument Type: Dissector
Length: 33cm 42cm
Curvature: Straight
Tip: Oval
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Breast Dissectors Angulated Blades Malleable

SKU: SSPS 315

Instrument Type: Dissector
Length: 33cm and 42cm
Size: 25mm Wide
Curvature: Angulated
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Breast Dissectors Angulated Blades Rigid

SKU: SSPS 316

Instrument Type: Dissector
Length: 33cm and 41cm
Size: 25mm Wide
Curvature: Angulated
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Breast Dissectors Angulated Rigid Shaft Insulated

SKU: SSPS 317

Instrument Type: Dissector
Length: 33cm and 41cm
Size: 25mm Wide
Curvature: Angulated
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Breast Elevator

SKU: SSPS 318

Instrument Type: Breast Elevator
Length: 28.00cm
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Iconoclast Breast Spreader

SKU: SSPS 319

Instrument Type: Breast Spreader
Size: Two Sizes Small 220mm and Large 295mm
Other Info: Blunt Straight Blade
Material: Japanese Stainless Steel
Reusable
Non-Sterile



McCOLLUM-DINGMAN Breast Dissector

SKU: SSPS 320

Instrument Type: Dissector

Length: 31cm and 41cm

Material: Japanese Stainless Steel

Reusable

Non-Sterile



Penn Breast Dissector

SKU: SSPS 321

Instrument Type: Dissector

Size: 234mm Long Width 25mm

Curvature: Straight

Material: Japanese Stainless Steel

Reusable

Non-Sterile



REYNOLDS Transaxillary Breast Dissector 40 cm

SKU: SSPS 322

Instrument Type: Dissector

Length: 31cm and 41cm

Material: Japanese Stainless Steel

Reusable

Non-Sterile



SILVERSTEIN Breast Dissector Blunt 355 mm

SKU: SSPS 323

Instrument Type: Dissector
Length: 35.5cm
Tip Configuration: Blunt
Curvature: Curved
Material: Japanese Stainless Steel
Reusable
Non-Sterile



SOLZ Breast Dissector

SKU: SSPS 324

Instrument Type: Dissector
Size: Length 34cm or 37cm
Material: Japanese Stainless Steel
Reusable
Non-Sterile




ACKERMAN Areola Marker

SKU: SSPS 325

For vertical mammoplasty procedures
Used to simplify the pattern marking as part of preparation for the procedure.

The Ackerman marker is available in 2 sizes 16cm circumference for large breast reductions and 14cm circumference standard size. Made with Japanese premium steel.



Castroviejo Surgical Caliper

SKU: SSPS 326

Instrument Type: Surgical caliper
Length: 8.5cm
Increments: 0.5mm
Maximum Measurement: 20mm
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Freeman Areola Cutter-Marker Standard

SKU: SSPS 327

Instrument Type: Areola Markers Single ended or Double ended
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Freeman Areola Cutter-Marker

SKU: SSPS 328

For accurate and precise measure and marking of the areola during breast reconstruction or augmentation procedures. Precise centrally located aperture for easy centering.

Can be used with or without methylene blue

Features

Ideal for accurate marking of the areola

Semi-sharp edge helps to lightly scribe the skin



Galvao Nipple Marker

SKU: SSPS 329

Instrument Type: Nipple Marker

Size: 40mm or 50mm

Material: Japanese Stainless Steel

Reusable

Non-Sterile



Solz Areola Marker

SKU: SSPS 330

Areola marker used to measure and mark the nipple-areola complex prior to making incisions in breast augmentation procedures and to make sure the NAC is centered appropriately. It has a standard center hole and is also fenestrated to provide additional vertical and horizontal markings.

Variety of inner diameter sizes, from 36mm to 45mm.

Tip: Blunt

Instrument Type: Areola Marker

Curvature: Round

Material: Japanese Stainless Steel

Reusable

Non-Sterile, Latex-Free



Halpern Breast Marker

SKU: SSPS 331

Marker for repositioning of the nipple with screw adjustment that fixes the angle to facilitate marking the contralateral breast
Useful tool in breast augmentation.
Allows for increased accuracy in marking the areola and the surrounding breast tissue.
The marker ranges from 38- 42mm diameter 60mm long shanks



McKissock Keyhole Breast Reduction Markers

SKU: SSPS 332


Surgical Instrument that is designed for breast reduction procedures and delineating areas of breast tissue to be excised
It accurately marks amount of tissue for removal.
Arms are graduated in centimeters for precise measurements.
Used together with Freeman Marker



Mammoplasty Caliper

SKU: SSPS 333

Instrument Type: Breast caliper/marker
Length: 28.6cm
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Breast Retractor Fiber-Optic Ferreira with Suction

SKU: SSPS 334

Instrument Type: Breast Retractor
Length of handle: Standard 13cm
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Freeman Flap Retractor

SKU: SSPS 335

Instrument Type: Flap - Breast Retractor
Size: Blade width 12.7cm and Blade Width 20.5cm
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Freeman Flap Retractor Fibre Optic

SKU: SSPS 336

Useful tool in breast reconstruction or augmentation procedures.
Universal adaptor that is compatible with standard light sources.
Ergonomically designed to minimizing user fatigue
The deep-set blade offers extra leverage during the operation
when exposure of deeper breast tissues is required Features
Universal adaptor
Ergonomically designed
Material: Japanese Stainless Steel



Marx Breast Retractor

SKU: SSPS 337

Instrument Type: Breast Retractor
 Length of handle: 27cm
 Size: 30x 125 mm and 30 x 140mm
 Material: Japanese Stainless Steel
 Reusable
 Non-Sterile



Tebbetts Retractor Fiber Optic Guide and Suction Blade

SKU: SSPS 338

Instrument Type: Breast Retractor
 Length: Tubes 16mm-40mm Blades 8cm -19cm
 Tip: Serrated
 Material: Japanese Stainless Steel
 Reusable
 Non-Sterile



Breast Retractor Tebbetts Double-Ended

SKU: SSPS 339

Double-ended breast retractor with serrated blades
 Ergonomically designed handle that offer comfort when using the retractor for extended periods.
 The smaller, shorter blade provides excellent exposure for initial access and dissection.
 The longer blade provides expanded surgical access.



Freeman Rake Retractor

SKU: SSPS 340

Four prong retractor designed for rhytidectomy and mammoplasty procedures

This retractor features two recessed inner prongs to provide more even retraction of the dissected issue.

Size: 18cm long and 22mm wide



Biggs Facu Life Retractor - Fiber Optic

SKU: SSPS 341

Instrument Type: Biggs Face Lift Retractor - Fiber Optic

Size: Length 229mm 5.5cm wide

Material: Japanese Stainless Steel

Steel Grade: Premium OR-Grade Precision Crafted from Quality Japanese Stainless Steel

Reusable

Non-Sterile



Arcus Marginalis Dissector, 24 cm S-Shaped, Curved

SKU: SSPS 342

Instrument Type: Arcus Marginalis Dissector

Size: Length 24cm 8mm Wide

Material: Japanese Stainless Steel

Steel Grade: Premium OR-Grade Precision Crafted from

Reusable

Non-Sterile



Daniel Endoscopic Forehead Nerve Dissector

SKU: SSPS 343

Instrument Type: Endoscopic Forehead Nerve Dissector
 Other Info Available curved right-left or half curved or straight
 LSize: 23.5cm
 Material: Japanese Stainless Steel
 Reusable
 Non-Sterile



D'Assumpcao Facelife Marker

SKU: SSPS 344

Instrument Type: Facelift Marker
 Size: Length - Overall 15.8cm Tip /Jaw 6mm
 Material: Japanese Stainless Steel
 Reusable
 Non-Sterile



Tebbetts Spatula Semi-Malleable

SKU: SSPS 345

Semi-malleable breast spatula
 Designed to facilitate multiple functions.
 The inside of the blade features atraumatic teeth which will not harm implant shells.
 Size: Length 32cm 1 x 6 inch blade



Breast Spatula With Handle

SKU: SSPS 346

Breast Spatula with Semi Malleable Blade
38mm Wide x 185mm Length Overall Length 315mm



Aston Face Lift SuperCut Scissors 16.5 cm

SKU: SSPS 347

Instrument Type: Face Lift Scissors
Length: 16.5cm
Other info: Semi-sharp outer edges. Curved. Serrated.
Working Surface Style Right and Right
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Castanares Face Lift Scissors

SKU: SSPS 348

Instrument Type: Face Lift Scissors
Length: 15.9cm
Other info: Serrated.
Curvation: Straight
Working Surface Style Right and Right
Material: Japanese Stainless Steel
Reusable
Non-Sterile



Gorney Scissors Curved SuperCut Scissors 20.5 cm

SKU: SSPS 349

Instrument Type: Face Lift Scissors
Length: 20.5cm
Material: Japanese Stainless Steel
Reusable
Non-Sterile



PEAK JOSEPH Face Lift Scissors

SKU: SSPS 350

Instrument Type: Face Lift Scissors

Length: 14.5cm


Other info: Serrated. Sharp/Sharp or Blunt/Blunt tips

Curvation: Straight or Curved

Material: Japanese Stainless Steel

Reusable

Non-Sterile





SUBLIME SURGICAL LIMITED
45 Waters Edge,
Canterbury CT1 1WX
United Kingdom
Tel: +447427677804
Web: www.sublimesurgical.com
E-mail: info@sublimesurgical.com