

STERILIZATION TRAYS

High Quality Instruments





General Sterilization Trays

01-11



Endoscopic Baskets

12-13



Universal Sterilization Boxes

14-18



PRF Box & Instruments

19-21




Dental Cassette Trays

22-37



Implantology / Bur holders

38-39



Orthodontic Steel Mirrors

40



Hollow Ware

41




Tray Accessories


42-44

General Sterilization Trays





Material: 304 (ASTM)
Steel sheet : 1.0mm



Perforation : Square pattern
Lid & bottom : 6.4 x 6.4mm
Side : 4.8 x 4.8mm



Easy locking System



Pull up & Stackable
Handles



Chamfer Corner



Round Corner

Squar pattern Chamfer corner

Art. No.	L	W	H
SNT-2020	540	240	50
SNT-2024	540	240	70
SNT-2028	540	240	90
SNT-2030	480	240	50
SNT-2034	480	240	70
SNT-2038	480	240	90
SNT-2040	380	250	50
SNT-2034	380	250	70
SNT-2038	380	250	90
SNT-2050	250	240	50
SNT-2054	250	240	70
SNT-2058	250	240	90
SNT-2060	250	120	50
SNT-2064	250	120	70
SNT-2068	250	120	90

Square pattern Round corner

Art. No.	L	W	H
SNS-2020	540	240	50
SNS-2024	540	240	70
SNS-2028	540	240	90
SNS-2030	480	240	50
SNS-2034	480	240	70
SNS-2038	480	240	90
SNS-2040	380	250	50
SNS-2034	380	250	70
SNS-2038	380	250	90
SNS-2050	250	240	50
SNS-2054	250	240	70
SNS-2058	250	240	90
SNS-2060	250	120	50
SNS-2064	250	120	70
SNS-2068	250	120	90


LIDs (optional)




Art. No.	L	W
SNT-2021	540	240
SNT-2031	480	240
SNT-2041	380	250
SNT-2051	250	240
SNT-2061	250	120

Art. No.	L	W
SNS-2021	540	240
SNS-2031	480	240
SNS-2041	380	250
SNS-2051	250	240
SNS-2061	250	120

*More sizes available on demand




Material: 304 (ASTM)
Steel sheet : 1.0mm



Perforation : Round pattern
8.0mm & 4.8mm dia.



Easy locking System



Pull up & Stackable
Handles



Chamfer Corner



Round Corner

Round pattern Chamfer corner

Art. No.	L	W	H
SNR-2020	540	240	50
SNR-2024	540	240	70
SNR-2028	540	240	90

SNR-2030	480	240	50
SNR-2034	480	240	70
SNR-2038	480	240	90

SNR-2040	380	250	50
SNR-2034	380	250	70
SNR-2038	380	250	90

SNR-2050	250	240	50
SNR-2054	250	240	70
SNR-2058	250	240	90

SNR-2060	250	120	50
SNR-2064	250	120	70
SNR-2068	250	120	90

Round pattern Round corner

Art. No.	L	W	H
SNU-2020	540	240	50
SNU-2024	540	240	70
SNU-2028	540	240	90

SNU-2030	480	240	50
SNU-2034	480	240	70
SNU-2038	480	240	90

SNU-2040	380	250	50
SNU-2034	380	250	70
SNU-2038	380	250	90

SNU-2050	250	240	50
SNU-2054	250	240	70
SNU-2058	250	240	90

SNU-2060	250	120	50
SNU-2064	250	120	70
SNU-2068	250	120	90


LIDs (optional)



Art. No.	L	W
SNR-2021	540	240
SNR-2031	480	240
SNR-2041	380	250
SNR-2051	250	240
SNR-2061	250	120

Art. No.	L	W
SNU-2021	540	240
SNU-2031	480	240
SNU-2041	380	250
SNU-2051	250	240
SNU-2061	250	120


*More sizes available on demand



Material: 304 (ASTM)




Round Corners
Perforation : 2.5mm 22(swg)



Mesh Type : Spotted
5.5 x 5.5 x 1.2mm



Easy to lock



Pull up & Stackable
Handles



Art. No.	L	W	H
SRCM-1050	540	255	50
SRCM-1070	540	255	70
SRCM-1100	540	255	100

SRCM-2050	485	255	50
SRCM-2050	485	255	60
SRCM-2070	485	255	70
SRCM-2100	485	255	100
SRCM-6003	480	300	64

SRCM-3050	405	255	50
SRCM-3070S	405	255	70
SRCM-3100	405	255	100

Art. No.	L	W	H
SRCM-4040	255	245	40
SRCM-4050	255	245	50
SRCM-4060	255	245	60


SRCM-4070	255	245	70
SRCM-4100	255	245	100
SRCM-5050	255	120	50
SRCM-4060	255	120	60
SRCM-5070	255	120	70

SRCM-5100	255	120	100
SRCM-6001	300	120	64


LIDs (optional)




Art. No.	L	W
SWD-9005	255	120
SWD-9004	255	245
SWM-9008	250	240
SWD-9001	540	255
SWD-9002	485	255
SWD-9003	405	255
SWM-9006	480	370
SWM-9007	480	250



Material: 304 (ASTM)




Chamfer Corners
Perforation : 2.5mm x 22(swg)



Mesh Type : Spotted
5.5 x 5.5 x 1.2mm



Easy to lock



Pull up & Stackable
Handles



Art. No.	L	W	H
SCCM-1050	540	255	50
SCCM-1070	540	255	70
SCCM-1100	540	255	100


SCCM-2050	485	255	50
SCCM-2070	485	255	70
SCCM-2100	485	255	100
SCCM-6005	480	300	60
SCCM-6010	480	250	60
SCCM-6012	480	120	60

Art. No.	L	W	H
SCCM-3050	405	255	50
SCCM-3070	405	255	70
SCCM-3100	405	255	100

SCCM-4050	255	245	50
SCCM-4070	255	245	70
SCCM-4100	255	245	100
SCCM-6001	300	240	60
SCCM-6020	250	240	60

SCCM-5050	255	120	50
SCCM-6025	250	120	60
SCCM-5070	255	120	70

LIDs (optional)




Art. No.	L	W
SLD-9005	255	120
SLD-9004	255	245
SLM-9008	250	240
SLM-6001	300	240
SLD-9001	540	255
SLD-9002	485	255
SLD-9003	405	255
SLM-9006	480	370
SLM-9007	480	250




Material: 304 (ASTM)



Single Frame
Round Corners



Mesh Type : Woven
5.5 x 5.5 x 1.0mm



Easy to lock



Pull up & Stackable
Handles

Deep drawn, stacking rim and electrolytically polished trays are equipped with collapsible grips



Art. No.	L	W	H
SDF-1050	540	255	50
SDF-1070	540	255	70
SDF-1100	540	255	100
SDF-2050	485	255	50
SDF-2060	485	255	60
SDF-2070	485	255	70
SDF-2100	485	255	100
SDF-6003	480	300	65
SDF-3050	405	255	50
SDF-3070	405	255	70
SDF-3100	405	255	100
SDF-8501	540	250	50
SDF-9000	480	250	50
SDF-9010	480	250	70
SDF-9020	480	250	100
SDF-9003	480	300	50
SDF-9030	480	250	60

Art. No.	L	W	H
SDF-4050	255	245	50
SDF-4060	255	245	60
SDF-4070	255	245	70
SDF-4100	255	245	100
SDF-6002	300	240	60
SDF-5050	255	120	50
SDF-5060	255	120	60
SDF-5070	255	120	70
SDF-5100	255	120	100
SDF-6001	300	120	64
SDF-9016	400	250	50
SDF-9106	400	250	70
SDF-9005	400	300	50
SDF-9004	400	370	50
SDF-9017	450	220	50
SDF-9100	250	240	50
SDF-9130	250	240	60


Art. No.	L	W	H
SDF-8500	500	150	100
SDF-8950	850	200	100
SDF-9102	355	100	90
SDF-9200	440	110	45
SDF-9210	440	110	65
SDF-9002	355	100	90
SDF-8900	410	260	110
SDF-8000	350	160	50
SDF-8005	260	190	50
SDF-8050	380	250	50
SDF-8100	380	190	85
SDF-9110	250	240	70
SDF-9107	250	190	50
SDF-9103	300	250	50
SDF-9600	650	650	30

LIDs (optional)




Art. No.	L	W
SWL-9001	540	255
SWL-9002	485	255
SWL-9006	480	370
SWL-9007	480	250


Art. No.	L	W
SWL-9003	405	255
SWL-9013	380	250
SWL-9005	255	120
SWL-9004	255	245
SWL-9008	250	240




Material: 304 (ASTM)




Double Frame protection
Round Corners



Mesh Type : Woven
5.5 x 5.5 x 1.0mm



Easy to lock



Pull up & Stackable
Handles

Deep drawn, stacking rim and electrolytically polished trays are equipped with collapsible grips



Art. No.	L	W	H
SWD-1050	540	255	50
SWD-1070	540	255	70
SWD-1100	540	255	100
SWD-2050	485	255	50
SWD-2060	485	255	60
SWD-2070	485	255	70
SWD-2100	485	255	100
SWD-6003	480	300	65
SWD-3050	405	255	50
SWD-3070	405	255	70
SWD-3100	405	255	100
SWM-8501	540	250	50
SWM-9000	480	250	50
SWM-9010	480	250	70
SWM-9020	480	250	100
SWM-9003	480	300	50
SWM-9030	480	250	60

Art. No.	L	W	H
SWD-4050	255	245	50
SWD-4060	255	245	60
SWD-4070	255	245	70
SWD-4100	255	245	100
SWD-6002	300	240	60
SWD-5050	255	120	50
SWD-5060	255	120	60
SWD-5070	255	120	70
SWD-5100	255	120	100
SWD-6001	300	120	64
SWM-9016	400	250	50
SWM-9106	400	250	70
SWM-9005	400	300	50
SWM-9004	400	370	50
SWM-9017	450	220	50
SWM-9100	250	240	50
SWM-9130	250	240	60


Art. No.	L	W	H
SWD-8500	500	150	100
SWD-8950	850	200	100
SWD-9102	355	100	90
SWD-9200	440	110	45
SWD-9210	440	110	65
SWM-9002	355	100	90
SWM-8900	410	260	110
SWM-8000	350	160	50
SWM-8005	260	190	50
SWM-8050	380	250	50
SWM-8100	380	190	85
SWM-9110	250	240	70
SWM-9107	250	190	50
SWM-9103	300	250	50
SWM-9600	650	650	30

LIDs (optional)



Art. No.	L	W
SWL-9001	540	255
SWL-9002	485	255
SWL-9006	480	370
SWL-9007	480	250


Art. No.	L	W
SWL-9003	405	255
SWL-9013	380	250
SWL-9005	255	120
SWL-9004	255	245
SWL-9008	250	240



Mesh Type : Woven
2 x 2 x 0.71mm(22-swg)
Material: 304 (ASTM)



Single Frame
Round Corners



Mesh Type : Woven
2 x 2 x 0.7mm




Easy to lock



Pull up & Stackable
Handles


Deep drawn and electrolytically polished trays



Art. No.	L	W	H	Lid
SMM-5500	40	40	20	Yes
SMM-5600	80	40	20	Yes
SMM-5650	160	40	20	Yes
SMM-5655	200	40	40	Yes
SMM-5656	200	60	40	Yes
SMM-5700	70	70	30	Yes
SMM-5800	80	80	40	Yes
SMM-5900	90	90	50	Yes
SMM-5906	80	40	40	Yes
SMM-5950	110	110	50	Yes
SMM-6010	150	110	20	Yes
SMM-6020	180	130	20	No
SMM-6030	210	170	20	No
SMM-6035	230	150	50	No
SMM-9101	250	240	50	No
SMM-9104	300	250	50	No
SMM-7000	200	110	40	No
SMM-7100	225	110	50	No
SMM-7105	235	120	60	No
SMM-7120	220	150	30	No
SMM-7200	270	135	40	No
SMM-7205	280	145	50	No

Art. No.	L	W	H	Lid
SMM-6040	250	200	20	No
SMM-6060	300	180	20	No
SMM-7055	155	155	80	No
SMM-9105	300	250	50	No
SMM-9301	400	200	50	No
SMM-6500	470	250	20	No
SMM-7300	400	250	50	No
SMM-8001	380	250	50	No
SMM-8901	410	260	110	No
SMM-6900	500	400	50	No
SMM-8051	380	250	50	No
SMM-9001	480	250	50	No
SMM-9011	480	250	70	No
SMM-9023	540	300	60	No
SMM-8100	380	190	85	No
SMM-9120	280	120	85	No
SMM-9012	480	200	50	No

Super Fine / Ultra Fine Mesh Trays



Deep drawn and electrolytically polished trays



Material: 304 (ASTM)




Round Corners



Easy to lock




Mesh Type : Woven
Super Fine Mesh
1 x 1 x 0.41mm(27-swg)



Mesh Type : Woven
Ultra Fine Mesh
0.7 x 0.7 x 0.25mm(32-swg)

Art. No.	L	W	H
SSM-5500	40	40	20
SSM-5600	80	40	20
SSM-5650	160	40	20
SSM-5655	200	40	40
SSM-5656	200	60	40
SSM-5700	70	70	30
SSM-5800	80	80	40
SSM-5900	90	90	50
SSM-5906	80	40	40
SSM-5950	110	110	50
SSM-6010	150	110	20
SSM-6020	180	130	20
SSM-6030	210	170	20
SSM-6035	230	150	50

Art. No.	L	W	H
SUM-5500	40	40	20
SUM-5600	80	40	20
SUM-5650	160	40	20
SUM-5655	200	40	40
SUM-5656	200	60	40
SUM-5700	70	70	30
SUM-5800	80	80	40
SUM-5900	90	90	50
SUM-5906	80	40	40
SUM-5950	110	110	50
SUM-6010	150	110	20
SUM-6020	180	130	20
SUM-6030	210	170	20
SUM-6035	230	150	50




Material : 304(ASTM)



Art. No.	L	W	H
STT-6518	650	178	135
STT-5718	570	178	135



Art. No.	L	W	H
STT-4733	470	330	70/82
STT-4725	470	250	70/82
STT-5330	530	300	70/82



Art. No.	L	W	H
STT-4750	600	300	250
STT-4752	600	300	300



Art. No.	L	W	H
STT-4740	400	300	100
STT-4742	400	300	150
STT-4746	500	300	100
STT-4748	500	300	200



Art. No.	L	W	H
SMB-4852	200	140	125
SMB-4854	180	110	100



Art. No.	L	W	H
SMB-4862	200	140	125
SMB-4864	180	110	100



Art. No.	L	W	H
SMB-4840	180	110	100
SMB-4842	200	140	125


Endoscope

To ensure the functionality and long life of delicate endoscope instruments is recommended to place these in our Endoscope baskets with silicone system for washing & sterilization and to eliminate metal to metal contact.


- ✓ Maximum Perforation are made with objective of reprocessing contaminated flexible endoscope.
- ✓ Perforated sides walls allows desired permeation of ethylene oxide assuring through its open during sterilization .
- ✓ Special Light Patterns, punching processed.
- ✓ Designing elegance manufacturing solution.
- ✓ Take control of your sterilization process and start getting things done.
- ✓ 304 grade Stainless steel provides aesthetic finish and corrosion protection.
- ✓ Includes a 4-years satisfaction guarantee
- ✓ Medical grade silicone stabilizes the instruments when the cassette is closed.
- ✓ Silicone Holders are offers 180 degrees and 50 cycles of protection.
- ✓ Burr free.
- ✓ Easy locking system, Button Lock/Clip Lock.
- ✓ Electrolytically polished.




Sheet: 304(ASTM)
1.0 mm



Perforation : Round pattern
8.0mm & 4.0mm dia.



Easy locking System




Art. No.	L	W	H
SNR-4011	290	80	45
SNR-4021	290	160	45


SNR-4031	450	80	45
SNR-4032	450	160	45

SNR-4052	550	150	65
SNR-4054	640	150	65


*More sizes available on demand



Perforation :
Square pattern 9.5 x 9.5mm



Easy locking System




Art. No.	L	W	H
SNT-4011	290	80	45
SNT-4021	290	160	45

SNT-4031	450	80	45
SNT-4032	450	160	45

SNT-4052	550	150	65
SNT-4054	640	150	65


*More sizes available on demand



EndoScopic Baskets



Mesh Type : Spotted
5.5 x 5.5 x 1.2 mm
Material: 304(ASTM)
Rods: 5.0mm,2.0mm



Mesh Type : Spotted
5.5 x 5.5 x 1.2mm




Easy to lock



Art. No.	L	W	H
SEB-4010	290	80	52
SEB-4030	460	80	52
SEB-4060	670	80	52
SEB-4070	700	80	52

SEB-4020	290	160	52
SEB-4120	310	150	77
SEB-4040	460	160	77
SEB-4050	640	150	77
SEB-4080	660	150	77
SEB-4110	460	160	52
SEB-4130	580	160	52


Mesh Type : Spotted
5.5 x 5.5 x 1.2 mm
Material: 304(ASTM)
Rods: 5.0mm,2.0mm



Art. No.	L	W	H
ECT-4010	290	80	52
ECT-4030	460	80	52
ECT-4020	290	160	52
ECT-4110	460	160	52

Universal Trays

Suitable for Surgical sets,
Orthopaedic and Ophthalmic Sets





Material : 304(ASTM)
Sheet: 1.0 / 1.2 mm



Art. No.	L	W	H	
SUT-2100	480	250	100	(3 layers)
SUT-3070	400	200	100	(3 layers)
SUT-3071	450	250	100	(3 layers)
SUT-6050	300	200	70	(2 layers)
SUT-6051	350	250	70	(2 layers)
SUT-4050	250	240	70	(2 layers)
SUT-4070	250	200	70	(2 layers)

Aluminum

Sheet : 1.5 mm



Art. No.	L	W	H	
SUA-2100	480	250	100	(3 layers)
SUA-3070	400	200	100	(3 layers)
SUA-3071	450	250	100	(3 layers)
SUA-6050	300	200	70	(2 layers)
SUA-4050	250	240	70	(2 layers)
SUA-4070	250	200	70	(2 layers)

*More sizes available on demand



Silicon Strips Fittings

Art. No.	L	W	H	
SDT-2121	500	250	100	(3 layers)
SDT-2123	400	200	100	(3 layers)
SDT-2126	350	250	60	(2 layers)
SDT-2129	350	250	100	(3 layers)
SDT-2131	250	240	60	(2 layers)

*More sizes available on demand



Silicon Strips Fittings

Art. No.	L	W	H	
SDT-3464	250	150	22	
SDT-3466	250	150	32	
SDT-3468	300	200	52	(2 layers)


*More sizes available on demand



Trays for Delicate Instruments



Art. No.	L	W	H
SNT-8003	275	175	30
SNT-8006	210	150	30



*More sizes available on demand




Art. No.	L	W	H
SMT-3464	250	150	22
SMT-3466	250	150	32
SMT-3468	300	200	52



Silicon mats is not included in price



Mesh Type : Spotted
5.5 x 5.5 x 1.2 mm
Material: 304(ASTM)
Sheet : 1.0 / 1.2 mm



Art. No.	L	W	H
SOT-2100	480	250	100
SOT-3070	400	250	70
SOT-6050	300	200	50
SOT-4050	250	240	50
SOT-4070	250	240	70



Art. No.	L	W	H
SMT-7040	300	200	40
SMT-7050	300	200	50
SMT-8040	400	250	40
SMT-8050	400	250	50



Art. No.	L	W	H
SMT-3456	300	200	40
SMT-3457	250	150	40



Needle Sterilization box



Art. No.	L	W	H
SYB-2901	150	90	16
SYB-2902	220	110	20



PRF System



PRF & GRF Box System Platelet Rich Fibrin Dental Implant Membrane Box


- ✓ Our Standardized and efficient preparation of PRF box Kit is designed to collect and transform PRF clots into membranes in sterile conditions and to conserve them in a clean and wet environment before use.
- ✓ The box also contains compression wells and maces to compress the PRF clots into dense PRF cylinders, easy to use for filling cavities.
- ✓ The serum exudate collected in the bottom of the box, can be used for a longer conservation of the membranes, and is ready to be mixed with a bone biomaterial for grafting.



SRF-9876
For Professional



SRF-9878
For Student



SRF-9880 09 - Holes
SRF-9881 12 - Holes



SRF-9883 12 - Holes



SRF-9885 12 - Holes



SRF-9888
68 x 60 x 20mm



SRF-9890
110 x 75 x 5 mm




SRF-9892
65 x 30 mm



SRF-9930 **SRF-9931** **SRF-9932** **SRF-9933** **SRF-9934**



SRF-9935 170 x 75 x 23 mm
cassette for soft brushing kit



Cassette for 7 - Instruments

Art. No.	L	W	H
SRF-9920	220	125	25



Cassette for PRF Instruments

Art. No.	L	W	H
SRF-9922	220	200	34

Cassette Trays




Cassette Trays:


- ✓ Our stainless-steel instrument cassettes offer proven reliability and quality performance with an unparalleled satisfaction guarantee.
- ✓ Our variety of Unique designs offers the years of trouble-free operation.
- ✓ Our Compact & durable features accommodates more instruments for greater efficiency in ultrasonic cleaning, autoclaving and storing.
- ✓ Apertures in lid and bottom of trays permit optimum flow of cleaning solution, ultrasonic energy and greater access for rinsing.
- ✓ Your instruments are completely protected throughout the use cycle with medical grade silicone instrument holders configured to the tops and bottoms of our trays.
- ✓ Various lengths, widths and heights of trays are available to accommodate different sizes and quantities of instruments within a particular Instruments.
- ✓ The durable, long-lasting, and soft medical grade silicone material provides maximum protection and extended use.
- ✓ Micro instruments cassettes & bur Holders are easy to organize in ultrasonic unit to remove bio-burden and to clean & protect instruments.
- ✓ Allows instruments & cassette/bur Holders organizer to dry thoroughly.
- ✓ Unique flow-through design allows for effective rinsing and cleaning.
- ✓ Short length hinge provides stable and secure attachment to cassette.
- ✓ Quick and easy button release provides a secure lock when closed.
- ✓ 304 grade Stainless steel provides aesthetic finish and corrosion protection.
- ✓ Includes a 4-years satisfaction guarantee
- ✓ Medical grade silicone stabilizes the instruments when the cassette is closed.
- ✓ Silicone Holders are offers 180 degrees and 50 cycles of protection.



Instruments Plus with adjustable driver & Detachable





Art. No.	L	W	H		
SST-2205	202	130	32	05 plus instruments	05 instruments + accessory area
SST-2208	202	140	32	08 plus instruments	08 instruments + accessory area
SST-2210	202	195	32	10 plus instruments	10 instruments + accessory area
SST-2212	267	165	32	12 plus instruments	12 instruments + accessory area
SST-2215	280	202	32	15 plus instruments	15 instruments + accessory area
SST-2216	280	202	32	16 plus instruments	16 instruments + accessory area
SST-2214	370	202	32	14 plus instruments	4 Elevator + 10 instruments + accessory area
SST-2220	370	202	32	20 plus instruments	20 instruments + accessory area



Art. No.	L	W	H		
SCT-2205	202	130	32	05 plus instruments	05 instruments + accessory area
SCT-2208	202	140	32	08 plus instruments	08 instruments + accessory area
SCT-2210	202	195	32	10 plus instruments	10 instruments + accessory area
SCT-2212	267	165	32	12 plus instruments	12 instruments + accessory area
SCT-2215	280	202	32	15 plus instruments	15 instruments + accessory area
SCT-2216	280	202	32	16 plus instruments	16 instruments + accessory area
SCT-2214	370	202	32	14 plus instruments	4 Elevator + 10 instruments + accessory area
SCT-2220	370	202	32	20 plus instruments	20 instruments + accessory area




Twin Cassette



Art. No.	L	W	H	
SCT-1910	200	85	44	10 instruments Detachable
SCT-1920	200	140	44	20 instruments Detachable

Also available with Hinges



Art. No.	L	W	H	
SST-1910	200	85	44	10 instruments Hinges
SST-1920	200	140	44	20 instruments Hinges



Also available with Detachable




Art. No.	L	W	H	
SCT-9810	200	182	34	For osteotome set
SCT-9813	200	200	34	For osteotome set with hammer




Bracket Lock



Art. No.	L	W	H	
SCT-1505	180	82	23	05 - Instruments
SCT-1507	180	105	23	07 - Instruments
SCT-1510	180	135	23	10 - Instruments
SCT-1520	265	180	23	20 - Instruments




Art. No.	L	W	H	
SKT-1505	190	70	23	05 - Instruments
SKT-1507	190	96	23	07 - Instruments
SKT-1510	190	132	23	10 - Instruments
SKT-1520	265	190	23	20 - Instruments




Art. No.	L	W	H	
SHT-1505	190	82	23	05 - Instruments
SHT-1507	190	105	23	07 - Instruments
SHT-1510	190	135	23	10 - Instruments
SHT-1520	265	190	23	20 - Instruments




Strip Lock




Art. No.	L	W	H	
SHC-1005	180	80	22	05 - Instruments
SHC-1007	180	105	22	07 - Instruments
SHC-1010	180	132	22	10 - Instruments
SHC-1020	260	180	22	20 - Instruments



Art. No.	L	W	H	
SJC-1005	180	80	22	05 - Instruments
SJC-1007	180	105	22	07 - Instruments
SJC-1010	180	130	22	10 - Instruments
SJC-1020	260	180	22	20 - Instruments



Art. No.	L	W	H	
SKC-1005	180	80	22	05 - Instruments
SKC-1007	180	105	22	07 - Instruments
SKC-1010	180	130	22	10 - Instruments
SKC-1020	260	180	22	20 - Instruments




Art. No.	L	W	H	
SAC-1005	180	80	22	05 - Instruments
SAC-1007	180	105	22	07 - Instruments
SAC-1010	180	130	22	10 - Instruments
SAC-1020	260	180	22	20 - Instruments




Strip Lock

Round bottom





Art. No.	L	W	H	
SCT-2105	185	90	23	05 - Instruments
SCT-2107	185	112	23	07 - Instruments
SCT-2110	185	145	23	10 - Instruments
SCT-2120	265	185	23	20 - Instruments




Art. No.	L	W	H	
SKT-2105	185	90	23	05 - Instruments
SKT-2107	185	112	23	07 - Instruments
SKT-2110	185	145	23	10 - Instruments
SKT-2120	265	185	23	20 - Instruments

Round bottom (Full Lid)





Art. No.	L	W	H	
SLT-2105	185	90	23	05 - Instruments
SLT-2107	185	112	23	07 - Instruments
SLT-2110	185	145	23	10 - Instruments
SLT-2120	265	185	23	20 - Instruments




Art. No.	L	W	H	
SKC-2105	185	90	23	05 - Instruments
SKC-2107	185	112	23	07 - Instruments
SKC-2110	185	145	23	10 - Instruments
SKC-2120	265	185	23	20 - Instruments



Square




Art. No.	L	W	H	
SCT-1205	180	70	22	05 - Instruments
SCT-1207	180	95	22	07 - Instruments
SCT-1210	180	130	22	10 - Instruments
SCT-1220	260	180	22	20 - Instruments




Art. No.	L	W	H	
SCT-9920	185	90	24	05 - Instruments
SCT-9921	185	112	24	07 - Instruments
SCT-9922	185	145	24	10 - Instruments
SCT-9923	265	185	24	20 - Instruments




Elongated



Art. No.	L	W	H	
SCT-2305	190	76	30	05 - Instruments
SCT-2307	190	102	30	07 - Instruments
SCT-2308	190	115	30	08 - Instruments
SCT-2310	190	142	30	10 - Instruments
SCT-2320	270	190	30	20 - Instruments
SCT-2216	280	202	30	20 - Instruments




Art. No.	L	W	H	
SCT-1605	200	80	35	05 - Instruments
SCT-1607	200	105	35	07 - Instruments
SCT-1610	200	145	35	10 - Instruments
SCT-1615	200	208	35	15 - Instruments
SCT-1620	280	200	35	20 - Instruments



Art. No.	L	W	H	
SYT-2705	222	95	22	05 - Instruments
SYT-2706	222	110	22	06 - Instruments
SYT-2707	222	126	22	07 - Instruments
SYT-2708	222	142	22	08 - Instruments
SYT-2710	222	174	22	10 - Instruments
SYT-2712	222	205	22	12 - Instruments
SYT-2714	268	202	30	12 - Instruments Plus
SYT-2720	370	230	32	20 - Instruments Plus




Elongated Single Lock



Art. No.	L	W	H	
SCT-2005	182	78	32	05 - Instruments
SCT-2007	182	102	32	07 - Instruments
SCT-2010	182	142	32	10 - Instruments
SCT-2020	270	182	32	20 - Instruments




Elongated (Double Lock)




Art. No.	L	W	H	
SHT-2003	203	50	32	03 - Instruments
SHT-2005	203	76	32	05 - Instruments
SHT-2007	203	102	32	07 - Instruments
SHT-2008	203	142	32	08 - Instruments + accessory area
SHT-2010	203	142	32	10 - Instruments
SHT-2012	265	165	32	12 - Instruments + accessory area
SCT-2216	280	202	30	20 - Instruments + accessory area



Elongated




Art. No.	L	W	H	
SCT-1403	180	55	23	03 - Instruments
SCT-1405	180	70	23	05 - Instruments
SCT-1407	180	95	23	07 - Instruments
SCT-1410	180	130	23	10 - Instruments
SCT-1420	260	180	23	20 - Instruments



Art. No.	L	W	H	
SYT-1403	190	55	23	03 - Instruments
SYT-1405	190	70	23	05 - Instruments
SYT-1407	190	95	23	07 - Instruments
SYT-1410	190	130	23	10 - Instruments
SYT-1420	260	190	23	20 - Instruments




UnionJack



Art. No.	L	W	H	
SCT-1805	185	70	22	05 - Instruments
SCT-1807	188	95	22	07 - Instruments
SCT-1810	185	130	22	10 - Instruments
SCT-1820	260	185	22	20 - Instruments


Cassette tray with steel fittings



Art. No.	L	W	H	
SYT-3105	185	80	23	05 - Instruments
SYT-3107	185	100	23	07 - Instruments
SYT-3110	185	135	23	10 - Instruments
SYT-3120	275	185	23	20 - Instruments



Double Hinge Cassette Trays



Art. No.	L	W	H	
SDT-2755	202	82	34	05 - Instruments
SDT-2756	202	90	34	06 - Instruments
SDT-2757	202	104	34	07 - Instruments
SDT-2760	202	142	34	10 - Instruments
SYT-2765	275	202	34	15 - Instruments & accessories area
SYT-2770	300	202	34	16 - Instruments & accessories area

Cassette tray Suitable for Photographic Mirrors




Art. No.	L	W	H
SMC-4202	140	80	15
SMC-4204	170	70	15
SMC-4206	200	80	15






Mesh Cassette



Mesh Type : Woven
5.5 x 5.5 x 1.0mm



Art. No.	L	W	H	
SMT-2505	195	80	32	05 - Instruments
SMT-2507	195	105	32	07 - Instruments
SMT-2510	195	140	32	10 - Instruments
SMT-2520	280	195	32	20 - Instruments



Art. No.	L	W	H	
SCT-9902	190	70	32	05 - Instruments
SCT-9903	190	102	32	07 - Instruments
SCT-9904	190	135	32	10 - Instruments



Cassettes for Orthodontic Pliers



Art. No.	L	W	H	
SPB-1224	185	58	78	04 - Orthodontic Pliers & accessories area
SPB-1226	185	85	78	06 - Orthodontic Pliers & accessories area
SPB-1230	185	130	78	10 - Orthodontic Pliers & accessories area



Art. No.	L	W	H	
SPB-1218	203	178	24	04 - Orthodontic Pliers & 4 Instruments
SPB-1220	203	240	24	06 - Orthodontic Pliers & 4 Instruments



Art. No.	L	W	H	
SPB-1210	187	150	70	for 10 Pliers



Art. No.	L	W	H	
SPB-1220	380	230	40	12-Pcs Extracting Pliers






Trays for Hand Piece & Small Parts




Art. No.	L	W	H
SMT-4004	80	40	20
SMT-4005	80	60	25
SMT-4006	100	60	25


Art. No.	L	W	H
SMT-4104	80	40	20
SMT-4105	80	60	25
SMT-4106	100	60	25



Art. No.	L	W	H
SMT-4012	150	40	30



Art. No.	L	W	H
SMT-4112	150	40	30



Art. No.	L	W	H
SMT-4113	150	40	30



Super Fine Mesh Trays



Ultra Fine Mesh Trays

Art. No.	L	W	H
SSM-5500	40	40	20
SSM-5700	70	70	30
SSM-5600	80	40	20
SSM-5906	80	40	40
SSM-5900	90	90	50
SSM-5650	160	40	20
SSM-5651	160	40	40

Art. No.	L	W	H
SUM-5500	40	40	20
SUM-5700	70	70	30
SUM-5600	80	40	20
SUM-5906	80	40	40
SUM-5900	90	90	50
SUM-5650	160	40	20
SUM-5651	160	40	40



Implant & Surgical bo



Art. No.	L	W	H
SIB-1010	100	80	48

Art. No.	
SEB-1009	1+8 Fittings



Art. No.	L	W	H
SIB-1012	180	120	52

Art. No.	L	W	H
SIB-1014	154	94	52

Art. No.	L	W	H
SIB-1016	180	120	52



Burr Holders

Art. No.	
SBH-1004	04 - Holes
SBH-1005	05 - Holes
SBH-1006	06 - Holes
SBH-1008	08 - Holes
SBH-1010	10 - Holes
SBH-1012	12 - Holes
SBH-1016	16 - Holes

SBH-1061	06 - Holes with Ratchet
-----------------	-------------------------



Burr Holders/Endo Box

Art. No.	
SBH-1106	07 - Holes
SBH-1107	09 - Holes
SBH-1112	12 - Holes
SBH-1015	15 - Holes




Burr Holders

Art. No.	
SBC-1004	04 - Holes
SBC-1005	05 - Holes
SBC-1006	06 - Holes
SBC-1008	08 - Holes
SBC-1010	10 - Holes
SBC-1012	12 - Holes
SBC-1016	16 - Holes




Burr Holders


Art. No.	
SBG-1006	06 - Holes
SBG-1008	08 - Holes
SBG-1012	12 - Holes
SBG-1016	16 - Holes




Orthodontics S. Steel Mirrors
Set of 7-
Front side reflection



IntraOral Photo Mirrors Stainless Steel



Photography Mirrors with handle Stainless Steel





Instruments Trays



HIT-8601	8" x 6" x 1-1/2"
HIT-1061	10" x 6" x 1-1/2"
HIT-1282	12" x 8" x 2"
HIT-1410	14" x 10" x 2"
HIT-1610	16" x 10" x 3"

available with LID, Perforation and Filter



Kidney Dish

HKD- 1006	Kidney Dish - 6"
HKD- 1008	Kidney Dish - 8"
HKD- 1010	Kidney Dish - 10"
HKD- 1012	Kidney Dish - 12"



Vomit Bowl



Lotion Bowl

HVB-1505	150 x 50mm
HLB-1205	205mm / 2 Ltr.
HLB-1150	150mm / 0.75 Ltr.
HLB-1125	125mm / 0.4 Ltr.



Gallipot

HGP-2002	2 oz (50ml)
HGP-2004	4 oz (100ml)
HGP-2006	6 oz (150ml)
HGP-2008	8 oz (200ml)



Dressing Container

HDC-1056	5" x 6-1/2"
HDC-1055	5" x 5"
HDC-1045	4" x 5"



Forceps Jar

HFJ-1025	2" x 5"
HFJ-1028	2" x 8"



Scaler Tray

HST-1901	190 x 100
-----------------	-----------

Sterilizing/Dressing Drum



	Ø	H
HSD-1510	150	100
HSD-1515	150	150
HSD-1919	190	190
HSD-1924	190	240
HSD-2419	240	190
HSD-2424	240	240
HSD-2439	240	390
HSD-2434	240	340
HSD-2915	290	150

Male Urinal



Female Urinal



HMU-1280	750ml
HFU -1280	750ml



Plastic Tray Tags
90 x 38

Identification Tags



SS Tag 50 x 30



SS Tag 90 x 30




Clip for cassettes

Adjustable Instruments Fixation Rack





Art. No.	L	W	H
STR-7002	230	220	40
STR-7004	200	190	40
STR-7006	220	80	40

Silicon Mats



Art. No.	L	W
SSM-1814	460	350 mm
SSM-1810	460	250 mm
SSM-1208	300	200 mm
SSM-1006	250	150 mm





Art. No.	
SCL-1002	05 - Instruments
SCL-1005	07 - Instruments
SCL-1008	10 - Instruments
SCL-1010	20 - Instruments

**Silicon
Color Codes**



A



B



C



D




E



F



G





Art. No.	
SCL-1012	05 - Instruments
SCL-1015	07 - Instruments
SCL-1018	10 - Instruments
SCL-1020	20 - Instruments



Art. No.	
SCL-1022	05 - Instruments
SCL-1025	07 - Instruments
SCL-1028	10 - Instruments
SCL-1030	20 - Instruments




Art. No.	
SCL-1032	05 - Instruments
SCL-1035	07 - Instruments
SCL-1038	10 - Instruments
SCL-1040	20 - Instruments



Art. No.	
SCL-1042	05 - Instruments
SCL-1045	07 - Instruments
SCL-1048	10 - Instruments
SCL-1050	20 - Instruments



Art. No.	
SCL-1052	05 - Instruments
SCL-1055	07 - Instruments
SCL-1058	10 - Instruments
SCL-1060	20 - Instruments



Art. No.
SCL-1103 10 - Instruments



Art. No.
SCL-1107 for Laryngoscope



Art. No.
SCL-1112 for Fixation Rack



Art. No.
SCL-1118 for Elevators



01



02



03



04



05



06

For Endo Scope Baskets

Art. No.	Art. No.
01- SCL-1122	04- SCL-1028
02- SCL-1024	05- SCL-1030
03- SCL-1026	06- SCL-1032

Silicon
Color Codes



A



B



C



D



E



F



G

01



02



03



04



06



07

For Implants & Burr Holder

Art. No.	Art. No.
01- SCL-1151	04- SCL-1058
02- SCL-1054	05- SCL-1060
03- SCL-1056	06- SCL-1070
	07- SCL-1072



SUBLIME SURGICAL LIMITED
45 Waters Edge,
Canterbury CT1 1WX
United Kingdom
Tel: +447427677804
Web: www.sublimesurgical.com
E-mail: info@sublimesurgical.com