



SERVING YOUR HEALTH WITH EXCELLENCE

SURGICAL INSTRUMENT

High Quality Instruments



Index

<i>Dissecting Scissors</i>	Page 01
<i>Metzenbaum Scissors, Mayo Scissors</i>	Page 02
<i>Iris Scissors, Doyen Scissors</i>	Page 03
<i>Utility Scissors, Enterotomy Scissors, Gradle Scissors, Halsted Mosquito Forceps</i>	Page 04
<i>Kelly Forceps, Rochester-Pean Forceps, Rochester-Ochsner Forceps</i>	Page 05
<i>Ergonomic Forceps</i>	Page 06
<i>Coverglass Forceps, Thumb Dressing</i>	Page 07
<i>Gerald & Locking Forceps, Tissue Forceps, Adson Forceps</i>	Page 08
<i>Adson Forceps, Splinter & Potts-Smith Forceps, Davis Forceps, Cushing Forceps</i>	Page 09
<i>Eye Dressing Forceps, Tweezers</i>	Page 10
<i>Saws, Dieffenbach Serrefine, Skull Breaker, Rachitome</i>	Page 11
<i>Sheas, Bone Cutting Forceps, Chisels</i>	Page 12
<i>Mallets, Needles Post Mortem, Needles Dissecting</i>	Page 13
<i>Probes, Scalpels & Handles, Knives, Brain Knife</i>	Page 14
<i>Spatula, Spinal Cord Remover, Graves Vaginal Speculum, Pederson Vaginal Speculum</i>	Page 15
<i>Weisman Graves Speculum, Weisman Pederson Vaginal Speculum</i>	Page 16
<i>Cusco Vaginal Speculum Center Screw</i>	Page 17

Dissecting Scissors

12-025 Size: 5 1/2" (14.0cm)

Straight, Sharp/Blunt Points

12-035 Size: 6 1/2" (16.5cm)

Straight, Sharp/Blunt Points

12-044 Size: 8" (20.3cm)

Straight, Sharp/Blunt Points



12-055 Size: 5 1/2" (14.0cm)

Straight, Sharp/Sharp Points

12-065 Size: 6 1/2" (16.5cm)

Straight, Sharp/Sharp Points

12-066 Size: 7" (17.8cm)

Straight, Sharp/Sharp Points

12-074 Size: 8" (20.3cm)

Straight, Sharp/Sharp Points



12-085 Size: 5 1/2" (14.0cm)

Straight, Slightly Blunt Tips

12-095 Size: 6 1/2" (16.5cm)

Straight, Slightly Blunt Tips

12-096 Size: 7" (17.8cm)

Straight, Slightly Blunt Tips

12-110 Size: 8" (20.3cm)

Straight, Slightly Blunt Tips



12-155 Size: 5 1/2" (14.0cm)

Curved, Sharp/Sharp Point

12-165 Size: 6 1/2" (16.5cm)

Straight, Slightly Blunt Tips



12-185 Size: 5 1/2" (14.0cm)

Curved, Blunt/Points

12-195 Size: 6 1/2" (16.5cm)

Curved, Blunt/Points

12-210 Size: 8" (20.3cm)

Curved, Blunt/Points



Metzenbaum Scissors

12-225 Size: 5 3/4" (14.6cm)

Straight, Narrow/Blunt Points

12-235 Size: 7" (17.8cm)

Straight, Narrow/Blunt Points



12-240 Size: 7" (17.8cm)

Curved, Narrow/Blunt Points

12-255 Size: 9" (22.9cm)

Curved, Narrow/Blunt Points

12-260 Size: 9" (22.9cm)

Curved, Narrow/Blunt Points



Mayo Scissors

12-305 Size: 5 1/2" (14.0cm)

Straight, Slightly Blunt Tips

12-310 Size: 5 1/2" (14.0cm)

Curved, Slightly Blunt Tips

12-315 Size: 6 3/4" (17.1cm)

Straight, Slightly Blunt Tips

12-320 Size: 6 3/4" (17.1cm)

Curved, Slightly Blunt Tips

12-325 Size: 9" (22.9cm)

Straight, Slightly Blunt Tips

12-330 Size: 9" (22.9cm)

Curved, Slightly Blunt Tips



Iris Scissors

12-410 Size: 4 1/2" (11.4cm)
Straight, Sharp Points

12-425 Size: 4 1/2" (11.4cm)
Curved, Sharp Points

12-185 Size: 4 1/2" (11.4cm)
Angled Over Edge, Sharp Points

12-186 Size: 4 1/2" (11.4cm)
Angled, Sharp/Blunt Points

12-190 Size: 4 1/2" (11.4cm)
Angled, Over Edge, Prob Points



Doyen Scissors

12-635 Size: 7" (17.8cm)
Straight, Blunt/Blunt Points



12-640 Size: 7" (17.8cm)
Curved, Blunt/Blunt Points



Utility Scissor

13-874 Size: 7 1/2" (19.1cm)
Serrated Shears, Plastic Handle Rings



Enterotomy Scissors

37-106 Size: 8 1/4" (20.3cm)
Straight, with one hook Blade



37-108 Size: 8 1/4" (20.3cm)
Straight, with one Probe Blade



Gradle Scissors

26-310 Size: 3 3/4" (9.4cm)
Slightly Curved, Sharp Points



Halsted Mosquito Forceps

14-055 Size: 4 7/8" (12.4cm)
Straight with Serrated Tips

14-060 Size: 4 7/8" (12.4cm)
Curved with Serrated Tips



Kelly Forceps

14-225 Size: 5 1/2" (14.0cm)
Straight Serrated Tips

14-230 Size: 5 1/2" (14.0cm)
Curved Serrated Tips




Rochester-Pean Forceps

14-405 Size: 5 1/2" (14.0cm)
Straight with Serrated Tips

14-425 Size: 7 1/4" (18.4cm)
Straight with Serrated Tips


14-445 Size: 9" (22.9cm)
Straight with Serrated Tips



Rochester-Ochsner Forceps

14-475 Size: 5 1/2" (14.0cm) (1x2 Teeth)
Straight with Serrated Tips

14-475 Size: 5 1/2" (14.0cm) (1x2 Teeth)
Curved with Serrated Tips



14-500 Size: 7 1/4" (18.4cm) (1x2 Teeth)
Curved with Serrated Tips

14-512 Size: 9" (22.9cm) (1x2 Teeth)
Straight with Serrated Tips

14-514 Size: 9" (22.9cm) (1x2 Teeth)
Curved with Serrated Tips



Ergonomic Forceps

15-500 Size: 7" (17.5cm)
Straight 3mm Smooth Tips

15-525 Size: 7" (17.5cm)
Straight 3mm Serrated Tips

15-545 Size: 7" (17.5cm)
Straight 3mm Tip 1x2 Teeth

15-550 Size: 7" (17.5cm)
Straight 1.5mm Smooth Tips

15-555 Size: 7" (17.5cm)
Straight 1.5mm Serrated Tips

15-560 Size: 7" (17.5cm)
Straight Smooth Needle Tips

15-575 Size: 7" (17.5cm)
Curved 1.5 Serrated Tips

15-580 Size: 7" (17.5cm)
Curved 1.5mm Smooth Tips


15-585 Size: 5 1/2" (14cm)
Straight 1.5mm Smooth Tips

15-590 Size: 5 1/2" (14cm) 30 Degree
Curved 1.5mm Smooth Tips

15-595 Size: 5 1/2" (14cm) 30 Degree
Curved 1.5mm Serrated Tips

15-600 Size: 5 1/2" (14cm)
Straight 1.5mm Serrated Tips

15-625 Size: 5 1/2" (14cm)
Straight 1.0mm Serrated Tips



Coverglass Forceps

11-102 Size: 4 1/4" (10.8cm) Angled Dual Enlarged Rectangular Tips	
26-286 Size: 4 1/4" (10.8cm) Curved Smooth Pointed Tips	

Thumb Dressing

17-030 Size: 4" (10.2cm) Straight Serrated Tips, Fine Grip 17-040 Size: 5" (12.7cm) Straight Serrated Tips, Fine Grip 17-045 Size: 5" (14.0cm) Straight Serrated Tips, Fine Grip 17-050 Size: 6" (15.2cm) Straight Serrated Tips, Fine Grip 17-065 Size: 12" (30.5cm) Straight Serrated Tips, Fine Grip	
17-047 Size: 6" (15.2cm) Straight, 2mm Smooth Tips, Broad Grip 17-049 Size: 7" (17.8cm) Straight, 2mm Smooth Tips, Broad Grip 17-051 Size: 6" (15.2cm) Straight, 2mm Serrated Tips, Broad Grip 17-052 Size: 7" (17.8cm) Straight, 2mm Serrated Tips, Broad Grip 17-055 Size: 8" (20.3cm) Straight, 2mm Serrated Tips, Broad Grip	 Fine Grip Broad Grip Smooth Tip Serrated Tip

Gerald & Locking Forceps

17-070 Size: 7" (17.8cm) Straight, Fine Grip, Serrated Tips	
17-075 Size: 7" (17.8cm) Curved, Fine Grip, Serrated Tips	
17-290 Size: 6" (15.2cm) Locking, Angled, Fine Grip, Serrated Tips	
17-291 Size: 6" (15.2cm) Locking, Angled, Fine Grip, Smooth Tips	

Tissue Forceps

17-130 Size: 4" (10.2cm) 1x2 Teeth, Tissue Forceps	
17-145 Size: 5 1/2" (14.0cm) 1x2 Teeth, Tissue Forceps	
17-152 Size: 7" (17.8cm) 1x2 Teeth, Tissue Forceps	
17-160 Size: 10" (25.4cm) 1x2 Teeth, Tissue Forceps	
17-350 Size: 6" (15.2cm) 3x4 Teeth, Tissue Forceps	
17-355 Size: 8" (20.3cm) 3x4 Teeth, Tissue Forceps	

Adson Forceps

17-555 Size: 4 3/4" (12.1cm) Serrated Tips	
17-556 Size: 6" (15.2cm) Serrated Tips	
17-557 Size: 7" (17.8cm) Serrated Tips	

Adson Forceps

17-560 Size: 4 3/4" (12.1cm)
1x2 Teeth



Splinter & Potts-Smith Forceps

18-030 Size: 6" (15.2cm)
Straight, Serrated Tips (A), Long Tip Length

18-031 Size: 6" (15.2cm)
Straight, Serrated Tips (A), Long Tip Length

23-356 Size: 7" (17.8cm)
Straight, 1.5mm Serrated Tips (B)

23-361 Size: 8 1/4" (20.9cm)
Straight, 1.5mm Serrated Tips

23-366 Size: 9 1/2" (24.1cm)
Straight, 1.5mm Serrated Tips



Davis Forceps

23-311 Size: 6" (15.2cm)
Straight, Serrated Spoon Type Tips

23-313 Size: 8" (20.3cm)
Straight, Serrated Spoon Type Tips

23-315 Size: 10" (25.4cm)
Straight, Serrated Spoon Type Tips



Cushing Forceps

24-181 Size: 7" (17.8cm)
Straight, Smooth Tips

24-182 Size: 7" (17.8cm)
Straight SP Type, Serrated Tips



Eye Dressing Forceps

26-275 Size: 4" (10.2cm)
Straight, Serrated Tips (A), Fine Grip

26-277 Size: 6" (15.2cm)
Straight, Smooth Tips, Broad Grip

26-278 Size: 7" (17.8cm)
Straight, Smooth Tips, Broad Grip

26-279 Size: 8" (20.3cm)
Straight, Smooth Tips (C), Broad Grip

26-280 Size: 4" (10.2cm)
Half Curved, Serrated Tips (D), Broad Grip



Tweezers

41-316
Overall Length: 3 1/2" (8.89cm) Tweezer Length: 3" (7.6cm)
Plastic Magnifier, Straight, Serrated Tips



41-317
Overall Length: 3 1/2" (8.89cm) Tweezer Length: 3" (7.6cm)
Plastic Magnifier, Curved, Serrated Tips



41-320
Overall Length: 5" (12.7cm) Tweezer Length: 4 1/2" (11.4cm)
Glass Magnifier, Straight, Smooth Tips, Broad Grip



Saws

13-420 Size: 10 1/2" (26.7cm) <i>Charriere</i>	
13-425 Size: 11 1/2" (29.2cm) <i>Satterlee</i>	
25-525 Size: 9" (22.9cm) <i>Langenbech Metacarpal Saw</i>	
13-426 Size: <i>Replacement Blade for Satterlee Saw</i>	

Dieffenbach Serrefine

19-725 Size: 1 3/4" (4.5cm) <i>Dieffenbach Serrefine Straight, Light Model</i>	
19-740 Size: 2 1/2" (6.4cm) <i>Dieffenbach Serrefine Curved, Light Model</i>	" For Pack of 6 " Visit Our Website 

Skull Breaker

37-635 Size: 5 1/2" (14.0cm) <i>Skull Breaker Hook End</i>	
37-636 Size: 5 1/2" (14.0cm) <i>Skull Breaker with Chisel End</i>	

Rachitome

37-660 Size: 8 3/4" (22.2cm) <i>Rachitome</i>	
---	---

Shears

23-610 Size: 8" (20.9cm)
Curved End, Plain Rib



23-635 Size: 8 3/4" (22.2cm)
Curved End, Gluck Rib



Bone Cutting Forceps


25-040 Size: 7 1/2" (19.1cm)
Straight, Ribbed Handles (A)

25-065 Size: 7 1/2" (19.1cm)
Curved, Ribbed Handles (B)

25-100 Size: 10" (25.4cm)
Curved, Ribbed Handles

25-090 Size: 7 1/2" (19.1cm)
Angular, Ribbed Handles (C)

25-075 Size: 10" (25.4cm)
Angular, Ribbed Handles



Chisels

25-605 Size: 6 1/4" (15.8cm) x 1/4" (6mm)
Army Pattern Chisels



25-610 Size: 6 1/4" (15.8cm) x 1/2" (12mm)
Army Pattern Chisels



25-615 Size: 6 1/4" (15.8cm) x 3/4" (18mm)
Army Pattern Chisels



Mallets

25-686 Size: 8 1/4" x 1" (21.6cm x 2.5cm) Mallet	
25-687 Size: 7 1/4" x 1" (17.8cm x 2.5cm) Mallet, Lead Filled	

Needles, Post Mortem

37-220 Size: 5" (12.7cm) Curved Post Mortem Needle	
37-222 Size: 4" (10.2cm) Curved Post Mortem Needle	
37-224 Size: 2 3/4" (7cm) Curved Post Mortem Needle	

" For Pack of 6 " Visit Our Website

Needles, Dissecting

37-505 Size: 5 1/4" (13.3cm) Teasing Needle Stainless Steel, Straight, Spade Tip, Black Plastic Handle	
37-510 Size: 5 1/4" (13.3cm) Teasing Needle, Stainless Steel, Straight, Pointed Tip, Black Plastic Handle	
37-550 Size: 5" (12.7cm) Dissecting Needle Stainless Steel, Straight, Pointed Tip, Metal Handle	
37-913 Size: 5 3/8" (13.6cm) Dissecting Needle Stainless Steel, Angled, Pointed Tip, Metal Handle	

Probes

27-510 Size: 5 1/2" (14.0cm) Double End Probe For Accurate Micro Autopsy Applications	
27-540 Size: 6" (15.2cm) Probe with Eyelet	
27-630 Size: 6" (15.2cm) Director, Grooved with Tongue Tie	 " For Pack of 6 " Visit Our Website

Scalpels & Handles

36-412 Size: 8 1/2" (21.6cm) Scalpel Handle #4L for #22 & #24 Surgical Blades	
36-419 Size: 5 1/2" (14.0cm) Scalpel Handle #4L for #22 & #22B Surgical Blades	

Knives

37-030 Cartilage Knife, Overall Length: 8" (20.3cm) Blade Length: 3 3/8" (8.6cm)	
37-035 Liston Post Mortem Knife, Overall Length: 12" (30.48cm) Blade Length: 3 3/8" (8.6cm)	
37-040 Liston Post Mortem Knife, Overall Length: 13.50" (34.29cm) Blade Length: 8" (20.3cm)	
37-060 Autopsy Knife, Overall Length: 8 1/4" (20.95cm) Blade Length: 3 1/4" (8.3cm)	

Brain Knife

23-315 Size: 10" (25.4cm) Straight, Serrated Spoon Type Tips	
--	--

Spatula

23-315 Size: 10" (25.4cm)
Straight, Serrated Spoon Type Tips



Spinal Cord Remover

23-315 Size: 10" (25.4cm)
Straight, Serrated Spoon Type Tips



Graves Vaginal Speculum

18-030 Size: 6" (15.2cm)
Straight, Serrated Tips (A), Long Tip Length

18-031 Size: 6" (15.2cm)
Straight, Serrated Tips (A), Long Tip Length

23-356 Size: 7" (17.8cm)
Straight, 1.5mm Serrated Tips (B)

23-361 Size: 8 1/4" (20.9cm)
Straight, 1.5mm Serrated Tips



Pederson Vaginal Speculum

18-030 Size: 6" (15.2cm)
Straight, Serrated Tips (A), Long Tip Length

18-031 Size: 6" (15.2cm)
Straight, Serrated Tips (A), Long Tip Length

23-356 Size: 7" (17.8cm)
Straight, 1.5mm Serrated Tips (B)

23-361 Size: 8 1/4" (20.9cm)
Straight, 1.5mm Serrated Tips



Weisman Graves Speculum

Right Side Open

50-15 Size: Medium
1 3/8" (3.5cm) x 4" (10.2cm)

50-20 Size: Large
1 1/2" (3.8cm) x 4 1/2" (11.4cm)

50-25 Size: Extra Large
1 1/2" (3.8cm) x 6" (15.2cm)

Left Side Open

50-30 Size: Medium
1 3/8" (3.5cm) x 4" (10.2cm)

50-35 Size: Large
1 1/2" (3.8cm) x 4 1/2" (11.4cm)

50-40 Size: Extra Large
1 1/2" (3.8cm) x 6" (15.2cm)



Weisman Pederson Vaginal Speculum

Left Side Open

60-55 Size: Medium
7/8" (2.2cm) x 4" (10.2cm)

60-60 Size: Large
1 2/5" (2.5cm) x 4 3/4" (12.1cm)

60-65 Size: Extra Large
1 1/2" (3.8cm) x 6" (15.2cm)

Right Side Open

60-70 Size: Medium
7/8" (2.2cm) x 4" (10.2cm)

60-75 Size: Large
1 2/5" (2.5cm) x 4 3/4" (12.1cm)

60-80 Size: Extra Large
1 1/2" (3.8cm) x 6" (15.2cm)




Cusco Vaginal Speculum Center Screw

70-70 Size: Extra Narrow
3/4" (1.9cm x 3 1/4" (8.3cm)

70-375 Size: Small
1 1/4" (3.2cm x 3 1/2" (8.9cm)

70-80 Size: Medium
1 3/8" (3.5cm x 4" (10.2cm)

70-85 Size: Large
1 3/8" (3.5cm x 4 1/2" (11.4cm)



Cusco Vaginal Speculum Side Screw

80-85 Size: Extra Narrow
3/4" (1.9cm x 3 1/4" (8.3cm)

80-90 Size: Small
1 1/4" (3.2cm x 3 1/2" (8.9cm)

80-95 Size: Medium
1 3/8" (3.5cm x 4" (10.2cm)

80-100 Size: Large
1 3/8" (3.5cm x 4 1/2" (11.4cm)





SUBLIME SURGICAL LIMITED
45 Waters Edge,
Canterbury CT1 1WX
United Kingdom
Tel: +447427677804
Web: www.sublimesurgical.com
E-mail: info@sublimesurgical.com